

HSW MULTI-MATIC®

Pistolet de vaccination pour bétail.

Repeater syringe for multiple dosing of live stock

25 ml / 50 ml

- Poignée durable en aluminium moulé sous pression
 - Une durabilité exceptionnelle, tout en étant léger
 - Cylindre en plastique medical durable
 - Avec raccord à aiguille métallique
 - biocompatible – autoclavable (134°C/273°F)
 - Résistant à tous les médicaments connus
 - Résistant aux UV et incassable
 - Dosage d'une précision exceptionnelle
 - tolérance < 1%
 - Réglage facile de la dose
 - via une molette de dosage de 0,5ml (25 ml) ou 1ml (50 ml)
 - Conception simplifiée
 - Peu de pieces de rechange nécessaires
 - Conception ergonomique pour une utilisation facile.
-
- Durable pressure-cast aluminium handle
 - exceptional durability, yet lightweight
 - Durable medical grade plastic barrel
 - With metal needle fitting
 - biocompatible – autoclavable (134°C/273°F)
 - chemically resistant to all known medication.
 - UV-resistant and break-proof
 - Exceptionally accurate dosing
 - tolerance < 1%
 - Easy dose setting
 - Via 0,5 ml (25 ml) or 1 ml (50 ml) dosage wheel
 - Streamlined design
 - Few spare parts required
 - Ergonomic design for ease of use




Version	Vol/ml	Ref.
Luer Lock	25 ml 50 ml	0840 1068



HENKE
SASS
WOLF

HSW MULTI-MATIC® 50ml

VETO CHIRURGICAL



Nr.	Description	Désignation	Ref.
2	Washer 50 ml	Rondelle 50 ml	1081
3	Piston+piston rod 50 ml	Tige de piston+piston 50 ml	1075
4	Locking sleeve cpl.	Tôle d'enclench. cpl.	1083
5	Spring for locking sleeve	Ressort tôle d'enclench.	1082
6	Knob black	Bouton noir	0186
7	Handle spring	Ressort de poignée	0189